FOREO



ISSA™ 2 sensitive set

THE WORLD'S 1ST SILICONE SONIC TOOTHBRUSH





FULL USER MANUAL

GETTING STARTED

Congratulations on taking the first step toward a whiter and brighter smile by acquiring this ISSA 2 *sensitive* set. Before you begin to enjoy all the benefits of sophisticated and stylish dental care in the comfort of your home, please take a few moments to look through this manual.

INTENDED USE: For in-home use, to promote overall oral hygiene by cleaning your teeth and gums.

Please **READ ALL INSTRUCTIONS BEFORE USE** and utilize this product only for its intended use as described in this manual.


MARNING: No modification of this equipment is allowed.

ISSA 2 sensitive set OVERVIEW

Offering the opportunity to experience the full ISSA 2 range, this exclusive set features the ISSA 2 *sensitive* 100% silicone brush head for the most gentle brushing experience and the ISSA 2 replaceable brush head for a more vigorous cleaning option. T-Sonic Pulsations in the toothbrush work together to remove plaque while gently massaging gums at 16 adjustable speeds and **2X the power**. With an exceptionally long battery life, a 1-hour USB charge powers up your device for **up to 365 days**, allowing for a whole-mouth clean feeling every day of the year!

GETTING TO KNOW YOUR ISSA 2 sensitive set

NOTE: The ISSA 2 *sensitive* may be locked when first unboxed – unlock it by holding down the + and – buttons together. The light will then flash to indicate that your device is unlocked.



FEATURES

ON/OFF Button:

- Press once to turn on, once to turn off
- After 3 minutes, your ISSA 2 sensitive will automatically turn off to prevent overuse


Intensity Controls:

- To lock and unlock the controls, hold down the + and buttons at the same time
- 16 intensities for maximum comfort: increase or decrease the intensity with the + or buttons
- Device will remember your preferences for your next use


Sensitive Brush Head - Soft silicone bristles for a more gentle brushing experience Hybrid Brush Head - Silicone and PBT polymer bristles are tough on plaque, gentle on gums Ultra-Hygienic Silicone Body - Nonporous to prevent bacteria build-up Indicator Light - Tells you when to switch to a new guadrant of the mouth & when to charge Charging Port - Seals the device making it 100% waterproof and usable in the shower USB Charger - Charge anytime, anywhere with the USB cable Travel Pouch - Protective sleeve keeps your ISSA 2 sensitive clean while travelling

HOW TO USE YOUR ISSA 2 sensitive set

Don't change your brushing routine, just your toothbrush!



1. Apply Toothpaste - Use the 2. Brush - Press the center button to suggested amount of your favorite toothpaste. For best results, we recommend using FOREO naturals Smile-Perfecting Toothpaste.



turn on your ISSA 2 sensitive and use the +/- buttons to adjust the intensity 4. Hybrid Brush Head (Included Atto your preferred setting. Brush your teeth using wide circular motions and light pressure for 2 minutes. ISSA 2 sensitive will prompt you to brush a new quadrant of the mouth every 30 seconds.




- 3. Rinse Turn off your device by pressing the center button. Rinse your mouth, and clean your ISSA 2 sensitive under running water. Repeat twice a day for a fresher, brighter smile!
 - tachment) For a more vigorous brushing experience, replace the sensitive brush head with the hybrid brush head and follow steps 1-3.

CAUTION: If you experience any discomfort when using your FOREO device, discontinue use immediately and consult a physician. If you find any blood on the brush surface, please clean the device with soap and water.

CLEANING YOUR ISSA 2 sensitive set

Clean the device after each use by washing the brush with warm water while running your fingers through the bristles. We recommend using FOREO Silicone Cleaning Spray on the body of your ISSA 2 *sensitive* and rinsing with warm water afterwards for optimal results. Allow device to air-dry. We do not recommend placing the ISSA 2 *sensitive* in a UV sanitation case.

NOTE: Never use cleaning products containing alcohol, petrol or acetone. Avoid using toothpastes with bleach on your device.



REPLACING YOUR ISSA 2 sensitive **BRUSH HEAD** PULL OFF, PUSH ON

All ISSA brush heads are compatible with the ISSA 2 range.

To replace a brush head, make sure your ISSA 2 *sensitive* is turned off, and lift the brush head straight up to detach. With the bristles of the replacement head facing forward, align its center hole with the metal dowel and push down to attach securely.

For maximum effectiveness, we recommend changing your ISSA 2 *sensitive* brush head every year. Each hybrid brush head lasts 6 months and Tongue Cleaner lasts up to 1 year.

You can order your replacement brush heads on www.foreo.com or by contacting the nearest authorized FOREO retailer.

CHARGING YOUR ISSA 2 sensitive

If the power of pulsations seems to diminish or the indicator light begins to blink, it's time to charge your device. Connect the USB charging cord to any outlet for 1 hour. When the indicator light stops blinking, your ISSA 2 *sensitive* is fully charged.

CAUTION: Before charging, make sure the charging port and charger are free of water and toothpaste residue. Do NOT use your FOREO device while charging, and do NOT charge for more than 24 hours. Only use the FOREO - provided USB charging cable.

IMPORTANT

FOR OPTIMUM SAFETY:

- Consult your physician prior to using the ISSA 2 sensitive if you have any medical concerns.
- If you are undergoing treatment for any oral condition, consult your dental professional prior to use.
- If you experience any discomfort when using your FOREO device, discontinue use immediately and consult a physician.
- Avoid leaving your device in direct sunlight and never expose it to extreme heat or boiling water.
- This device is not intended for children under the age of 5. Close supervision is necessary when this device is used by, on, or near children. Proper supervision and instructions for safe use should be given to children, as well as those with reduced physical and mental abilities.
- Children should not be allowed to play with the device.
- Discontinue use if this product appears damaged in any way (brush head, handle, charger). This product contains no serviceable parts.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Never use the device if it is damaged.

TROUBLESHOOTING

Precautions to be taken in the event of changes in the performance of the ISSA 2 sensitive:

The ISSA 2 sensitive is not activated when pressing the center button?

• Battery is empty. Recharge your ISSA 2 sensitive.

The ISSA 2 sensitive cannot be switched off and/or interface buttons do not respond?

• Microprocessor is temporarily malfunctioning. Connect the charging cable to restore function to normal.

The ISSA 2 sensitive's light does not blink when charger is connected?

- Battery is fully charged, ready for up to 365 days of brushing.
- Battery is drained and will require several minutes to acknowledge charging cable.
- Charging cable is not properly connected. Check electrical outlet/USB socket and charging port.

The ISSA 2 sensitive has fallen into water and the charging port is wet?

• Dry off and continue use. The ISSA 2 sensitive is completely waterproof and the charging port is fully sealed.

If you do not find the answer to your specific question in the troubleshooting section, or if you have any other questions regarding the device's operation, please visit the Customer Care section of www.foreo.com/support.

WARRANTY TERMS AND CONDITIONS

REGISTER WARRANTY

To activate your 2-Year Limited Warranty register the number provided on the FOREO Magnetic Scratch Card at www.foreo.com/support.

2-YEAR LIMITED WARRANTY

FOREO warrants this device for a period of TWO (2) YEARS after the original date of purchase against defects due to faulty workmanship or materials arising from Normal Use of the device. The warranty covers working parts that affect the function of the device. It does NOT cover cosmetic deterioration caused by fair wear and tear, or damage caused by accident, misuse or neglect. Any attempt to open or take apart the device (or its accessories) will void the warranty.

If you discover a defect and notify FOREO during the warranty period, FOREO will, at its discretion, replace the device free of charge. Claims under warranty must be supported by reasonable evidence that the date of the claim is within the warranty period. To validate your warranty, please keep your original purchase receipt together with these warranty conditions for the duration of the warranty period.

To claim your warranty, you must log in to your account at www.foreo.com and then select the option to make a warranty claim. Shipping costs are nonrefundable. This undertaking is in addition to your statutory rights as a consumer and does not affect those rights in any way.

*2-Year Limited Warranty does not apply to Replacement Brush Heads.

DISPOSAL INFORMATION

Disposal of old electronic equipment (applicable in the EU and other European countries with separate waste collection systems).

The crossed-out dustbin symbol indicates that this device should not be treated as household waste, but rather be brought to the appropriate collection point for recycling of electrical and electronic equipment. By ensuring this device is disposed of correctly, you will help prevent the potential negative consequences for the environment and human health which could be caused by inappropriate waste handling of the product. The recycling of materials will also help conserve natural resources.

For more information about the recycling of your device, please contact your local household waste disposal service or your place of purchase.



BATTERY REMOVAL:

NOTE: This process is not reversible. Opening the device will void its warranty. This action must only be undertaken when the device is ready to be disposed of.

Because this device contains a lithium-ion battery, the battery must be removed before disposal and should not be thrown away with household waste. To remove the battery, pull off the brush head, peel away the silicone cover, open the plastic shell and remove the battery to be disposed of in accordance with your local environmental regulations. Wear gloves during this process for your safety. Detailed visual instructions are provided below:



SPECIFICATIONS

MATERIALS:

COLORS: SIZE: WEIGHT: BATTERY: USER TIME: STANDBY: FREQUENCY: INTERFACE: PULSATIONS/MIN: Medical-grade Silicone + PBT Polymer Pearl Pink/Mint 195 × 32.4 × 115 mm 77g Li-Ion 1000mA 3.7V Up to 365 days 180 days 166 Hz 3 buttons Up to 11,000

ENVIRONMENTAL CONDITIONS FOR USE:

Temperature:	5° to 40° Celsius
Humidity:	40% to 80%
Atmospheric Pressure:	800 to 1,060 hPa

ENVIRONMENTAL CONDITIONS FOR STORAGE AND SHIPPING:

Temperature:	-10° to 50° Celsius
Humidity:	30% to 80%
Atmospheric Pressure:	500 to 1,060 hPa

DISCLAIMER:

Users of this device do so at their own risk. Neither FOREO nor its retailers assume any responsibility or liability for any injuries or damages, physical or otherwise, resulting, directly or indirectly, from the use of this device. Further, FOREO reserves the right to revise this publication and to make changes from time to time in the contents thereof without obligation to notify any person of such revision or changes.

Model may be changed for improvements without notice.

©2020 FOREO AB. ALL RIGHTS RESERVED.

EU DISTR.: FOREO AB, KARLAVÄGEN 41, 114 31 STOCKHOLM, SWEDEN. US DISTR.: FOREO INC., 1525 E PAMA LANE, LAS VEGAS, NV 89119, USA. AUS DISTR.: 6 PROHASKY STREET, PORT MELBOURNE, VIC, 3207 AUSTRALIA. DESIGNED AND DEVELOPED BY FOREO SWEDEN. MANUFACTURED FOR FOREO AB.

WWW.FOREO.COM