

FOREO




UFO™ 2

SWEDISH BEAUTY-TECH
POWER MASKS



ONLINE MANUAL

FULL USER MANUAL



GETTING STARTED

Congratulations on taking the first step toward out-of-this-world skincare by acquiring UFO™ 2. Before you begin to enjoy all the benefits of sophisticated power mask technology in the comfort of your home, please take a few moments to carefully read the instructions in this manual.


Please READ ALL INSTRUCTIONS BEFORE USE and utilize this product only for its intended use as described in this manual.

⚠ WARNING: NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED.

UFO™ 2 OVERVIEW

Re-engineered with ultra-fast heating and an even wider spectrum of LED light, UFO™ 2 is a revolutionary power mask device that gives you superior skin in seconds. Featuring advanced Hyper-Infusion Technology, UFO™ 2 combines Thermo-Therapy and Cryo-Therapy with T-Sonic™ pulsations to enhance the benefits of your mask. Further equipped with 8 different LED lights, UFO™ 2 helps you unmask a refined, radiant complexion in only 90 seconds!

GETTING TO KNOW YOUR UFO™ 2




UFO™ POWER ACTIVATED MASK
Soft microfiber masks exclusively formulated for use with UFO™.



USB CHARGING CABLE
Charge anytime, anywhere with USB cable.




INTERACTIVE APP
Automatically syncs pre-programmed power mask routines to device.



SLEEK & STYLISH STAND
Protects and displays your UFO™ 2, while helping it dry quickly.

FIRST THINGS FIRST...

- 1 To unlock and register your device for first time use, download the FOREO For You app.
- 2 Connect to the app via Bluetooth on your smartphone, and press the universal power button on UFO™ 2 to sync your device to the app.
- 3 In the app, scan the mask barcode and follow the instructions provided.



Download the FOREO app here



Remove the attachment ring from UFO™ 2.



Remove UFO™ Power Activated Mask from sachet.



Secure mask on device by clipping the attachment ring in place, with FOREO logo facing up.



Enjoy your power mask treatment!

HOW TO USE UFO™ POWER ACTIVATED MASKS

Start with clean, dry skin. We recommend cleansing with LUNA™ first for best results.

- 1 With a mask secured to your device, gently glide UFO™ 2 across your skin, distributing the mask essence evenly.
- 2 Using the UFO™ 2 device, massage the formula into your skin with a circular motion. Move the device in a clockwise circle around the face, starting on the chin, then moving up to the right cheek and the forehead. Continue down the left cheek to the chin and neck. Repeat this movement until the treatment ends.
- 3 Gently pat any remaining essence into your skin, or remove excess with a cotton pad. Apply a moisturizer of your choice, if desired.

*For the eye mask, gently glide around the entire eye contour, alternating between the left and right eye until the treatment ends.



FOREO SHEET MASKS

- 1 Gently unfold your FOREO sheet mask and place it on your face, adjusting around the eyes and mouth.
- 2 Using the UFO™ 2 device, massage the formula into your skin with a circular motion. Move the device in a clockwise circle around the face, starting on the chin, then moving up to the right cheek and the forehead. Continue down the left cheek to the chin. Repeat this movement until the treatment ends.
- 3 Leave the sheet mask on for another 10 minutes, or dispose of it as you prefer. Gently pat any remaining essence into your skin, or remove excess with a cotton pad. Apply a moisturizer of your choice, if desired.

You can now use your UFO™ 2 with most sheet masks, thanks to new Advanced Temperature Control!*

***NOTE:** Heat and cold affect different skincare products in different ways, and can compromise the ingredients in the formulas. For optimal safety, we strongly recommend researching the effects of temperature on any non-FOREO mask before enhancing it with UFO™ 2.

ADDING A NEW TREATMENT

UFO™ 2 can be programmed with up to 8 mask treatments! Make My Day and Call It a Night are pre-programmed to the first two indicator lights, and the third light is reserved for your customized treatment.

Any other mask treatment must be activated through the FOREO For You app for the first time. Then, the routine can be accessed through the app or by using the universal power button.



USING UFO™ 2 WITHOUT THE APP

- 1 Press the universal power button to turn on and activate the Make My Day treatment routine.
- 2 Press again to activate the Call It a Night treatment, and so on with each mask. (For example, if H2Overdose is on slot 4 in your app, press the universal power button 4 times to activate the H2Overdose treatment).
- 3 Your power mask treatment is ready to begin!

STILL WANT MORE?

If one UFO™ 2 treatment wasn't enough, press the universal power button again within 30 seconds after the end of the routine to start the same treatment over again! You can also personalize your power mask experience by creating your own treatment routine via the app, where you can select / adjust each technology to accompany your mask.

CLEANING UFO™ 2

Always clean UFO™ 2 thoroughly after use. Wash the device and plastic ring with soap and then rinse with warm water. Pat dry with a lint-free cloth or towel. After use, we recommend spraying the device with FOREO's Silicone Cleaning Spray and rinsing with warm water for optimal results, then placing the device on the stand to dry.

NOTE: Never use cleaning products containing alcohol, petrol or acetone, as they may irritate the skin and damage the silicone.

IMPORTANT FOR OPTIMUM SAFETY:

- If you have a skin condition or any medical concerns, please consult a dermatologist before use.
- Using UFO™ 2 should be comfortable – if you experience discomfort, discontinue use immediately and consult a physician.
- Exercise particular care when using UFO™ 2 around the under-eye regions and do not bring the device into contact with the eyelids or eyes themselves.
- At no time should the device be placed directly onto the skin/face with power on.
- For reasons of hygiene, we do not recommend sharing your UFO™ 2 with anyone else.
- Avoid leaving your UFO™ 2 in direct sunlight, and never expose it to extreme heat or boiling water.
- UFO™ 2 has a heated surface. Those sensitive to heat must use caution and care when using this device.
- FOREO has not evaluated and tested the safety and efficacy of UFO™ 2 devices in conjunction with all non-FOREO sheet masks and/or skincare products. Use this product with non-FOREO masks and skincare products at your own risk.
- Close supervision is necessary when this device is used by, on, or near children, as well as those with reduced physical and mental abilities.
- Discontinue use if this product appears damaged in any way. This product contains no serviceable parts.
- Use this device only for its intended use as described in this manual. If you do not find the answer to your specific question, or if you have any other questions regarding the device's operation, please visit www.foreo.com

TROUBLESHOOTING

Precautions to be taken in the event of changes in the performance of UFO™ 2.

If UFO™ 2 is not activated when pressing the universal power button:

- Battery is empty. Charge your UFO™ 2 via USB charging cable until the light emits a steady glow. A full 2.5 hour charge will last up to 40 treatments.

If UFO™ 2 cannot be switched off and/or interface buttons do not respond:

- Microprocessor is temporarily malfunctioning. Press and hold universal power button to restart the device.

If UFO™ 2 won't sync to the FOREO For You app:

- Switch your Bluetooth off and then on again to try reconnecting.
- Close the FOREO For You app and then reopen it to start the process over.
- Check to see if the app needs to be updated in your mobile device's app store.

WARRANTY TERMS & CONDITIONS

REGISTER WARRANTY

To activate your 2-Year Limited Warranty, register by downloading the FOREO For You app and following the instructions provided.

2-YEAR LIMITED WARRANTY

FOREO warrants this device for a period of TWO (2) YEARS after the original date of purchase against defects due to faulty workmanship or materials arising from Normal Use of the device. The warranty covers working parts that affect the function of the device. It does NOT cover cosmetic deterioration caused by fair wear and tear, or damage caused by accident, misuse or neglect. Any attempt to open or take apart the device (or its accessories) will void the warranty.

If you discover a defect and notify FOREO during the warranty period, FOREO will, at its discretion, replace the device free of charge. Claims under warranty must be supported by reasonable evidence that the date of the claim is within the warranty period. To validate your warranty, please keep your original purchase receipt together with these warranty conditions for the duration of the warranty period.

To claim your warranty, you must log in to your account at www.foreo.com and then select the option to make a warranty claim. Shipping costs are nonrefundable. This undertaking is in addition to your statutory rights as a consumer and does not affect those rights in any way.

DISPOSAL INFORMATION


Disposal of old electronic equipment (applicable in the EU and other European countries with separate waste collection systems).



The crossed-out dustbin symbol indicates that this device should not be treated as household waste, but rather be brought to the appropriate collection point for recycling of electrical and electronic equipment. By ensuring this device is disposed of correctly, you will help prevent the potential negative consequences for the environment and human health which could be caused by inappropriate waste handling of the product. The recycling of materials will also help conserve natural resources.

For more information about the recycling of your device, please contact your local household waste disposal service or your place of purchase.

BATTERY REMOVAL



SPECIFICATIONS

MATERIALS:	Body-safe silicone & PC + ABS, Aluminium Alloy	TEMPERATURE:	5 - 45 Degrees Celsius
COLOR:	Pearl Pink/Mint/Fuchsia/Black	USAGE:	Up to 40 uses per charge
SIZE:	72.5 mm x 31.4 mm	STANDBY:	180 days
WEIGHT:	146.2g	FREQUENCY:	180 Hz
BATTERY:	Li-Ion 1200mAh 3.7V	MAX NOISE LEVEL:	<50 dB
		INTERFACE:	1-button

Disclaimer: Users of this device do so at their own risk. Neither FOREO nor its retailers assume any responsibility or liability for any injuries or damages, physical or otherwise, resulting directly or indirectly from the use of this device. Furthermore, FOREO reserves the right to revise this publication and to make changes from time to time in the contents thereof without obligation to notify any person of such revision or changes.

CAUTION: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio / TV technician for help.

Model may be changed for improvements without notice.

FAQs

UFO™ 2 BASICS

1. WHAT SHOULD I DO AFTER I RECEIVE MY UFO™ 2 DEVICE?

Congratulations on discovering smarter skincare! First things first - download the FOREO For You app for free, to unlock and register your device. (For more information, please refer to the section below titled 'THE APP').

2. HOW DO I START MY FIRST TREATMENT?

First, clean and dry your face and neck carefully - we recommend cleaning with LUNA™ first for best results. Then attach a UFO™ Power Activated Mask to your UFO™ 2 device, or place a sheet mask on your face. Finally, activate your UFO™ 2 mask routine via the app, and follow the step-by-step instructions available. For more information, please refer to the How To Use guide above.

3. HOW DO I TURN MY UFO™ 2 DEVICE ON / OFF?

Press the universal power button on your UFO™ 2 device to turn it on. To turn your device off, press and hold the universal power button for 3 seconds. After completing a pre-programmed mask routine, the device will turn off automatically.

4. WHAT COMES WITH MY UFO™ 2 DEVICE?

1 UFO™ 2 Device, 1 USB Charging Cable, 1 Make My Day Mask, 1 User Manual, 1 Quick Start Guide, & 1 Stand.

THE UFO™ 2 DEVICE

5. HOW OFTEN CAN I USE MY UFO™ 2 DEVICE?

UFO™ 2 skincare routines have been developed for twice-a-day use, morning and evening. However, depending on your skin's needs and your personal preferences, you can use it more or less.

6. HOW LONG DOES EACH TREATMENT TAKE?

The Make My Day & Call It a Night mask treatments take just 90 seconds. All other treatment times will vary depending on the most effective combination of technologies.

7. IS EACH UFO™ 2 TREATMENT THE SAME?

No, each FOREO mask pairs with a unique treatment routine for optimal results.

8. IS MY UFO™ 2 DEVICE SAFE TO USE ON SENSITIVE SKIN?

UFO™ 2 devices are safe to use and suitable for all skin types. Furthermore, depending on your skin's specific needs, you can choose from a range of different FOREO Masks and customize your own treatment routine to perfectly cater to your skincare desires.

9. CAN I USE MY UFO™ 2 DEVICE AROUND THE EYE AREA?

Yes, UFO™ 2 devices are safe to be used under the eyes. We recommend gently gliding the device around the eye area, being careful to apply minimal pressure.

10. CAN I USE MY UFO™ 2 DEVICE AFTER HAVING PLASTIC SURGERY, BOTOX OR RADIESSE?

UFO™ 2 devices offer a gentle facial treatment, suitable for all skin types. However, you should consult your doctor before use if you have had any recent medical procedures.

11. CAN I USE MY FOREO LUNA™ WITH MY UFO™ 2 DEVICE?

Yes! UFO™ & LUNA™ devices work better together - LUNA™ facial cleansing brushes gently exfoliate dead skin cells, unclog pores and remove up to 99.5% of dirt and oil, making your skin the perfect canvas for a UFO™ 2 power mask treatment.

12. WHAT IS THE DIFFERENCE BETWEEN UFO™ 2 AND UFO™ mini 2?

Both UFO™ 2 and UFO™ mini 2 offer a sophisticated power mask treatment for you to enjoy in only 90 seconds. However, UFO™ 2 has an additional Cryo-Therapy feature, which instantly cools to shrink the appearance of pores and provide a firmer looking complexion. Also, while UFO™ mini 2 can only be used with selected FOREO masks, UFO™ 2 can be used with all FOREO masks and most non-FOREO sheet masks.*

13. WHAT ARE THE MAIN DIFFERENCES BETWEEN UFO™ DEVICES AND UFO™ 2 DEVICES?

Unlike UFO™ devices, UFO™ 2 devices have been re-engineered with ultra-fast Thermo-therapy and Cryo-therapy, to offer a more efficient and effective masking experience than ever before. UFO™ 2 also features advanced temperature control, allowing you to adjust the temperature of your Thermo- and Cryo-Therapy treatments to your personal preference. As a result, you can now use UFO™ 2 with most sheet masks.*

THE APP

14. WHY DO I NEED TO DOWNLOAD THE FOREO APP BEFORE USING MY UFO™ 2 DEVICE?

All app-enabled FOREO products need to be activated and unlocked via the FOREO For You app the first time they are turned on. This helps protect against fake/copied FOREO products, and makes it easier to register the warranty and protect your investment.

15. HOW DO I PAIR MY UFO™ 2 DEVICE WITH THE FOREO APP?

Download the FOREO For You app to your smartphone or tablet and turn on Bluetooth. Press the universal button on your device once to turn on. A blinking white light indicates your device is in pairing mode. Follow the instructions in the app to register and pair your UFO™ 2 device.

16. DOES MY UFO™ 2 DEVICE WORK WITHOUT THE APP?

Yes, you can manually use both devices and they do not require use of the FOREO For You app. However, to unlock your device for your first use, you must register your device through the FOREO For You app and select your preferences. The Make My Day and Call It a Night treatments are pre-programmed into UFO™ 2 and do not require use of the FOREO app. However, to unlock all remaining FOREO masks, you must either select your treatment or scan the barcode on the back of each FOREO mask sachet via the app to activate the treatments - the mask will then automatically sync with your device so it can be accessed manually next time. To customize your power mask routine, you must go directly through the app.

17. CAN I CUSTOMIZE MY POWER MASK TREATMENT ROUTINE?

Yes. While each UFO™ Power Activated Mask and FOREO Sheet Mask comes with a recommended pre-programmed treatment routine, you can personalize your power Mask experience by creating your own treatment routine via the FOREO For You app. On the app, you are able to choose your preferred LED light and T-Sonic™ pulsation level, and adjust Thermo-Therapy and Cryo-Therapy according to your preferences.

TROUBLESHOOTING & MAINTENANCE

18. WHAT DOES IT MEAN WHEN MY UFO™ 2 DEVICE IS CONTINUOUSLY BLINKING?

A blinking light around the device indicates that your UFO™ 2 device is in Bluetooth pairing mode or that the device needs to be charged. Follow the instructions in the app to sync your device.

19. HOW DO I CHARGE MY UFO™ 2 DEVICE?

UFO™ 2 devices are USB-rechargeable, with a two and half hour charge delivering up to 40 uses.

20. HOW DO I CLEAN MY UFO™ 2 DEVICE?

Wash the device and plastic ring with soap and then rinse with warm water. Pat dry with a lint-free cloth or towel. After use, we recommend spraying the device with FOREO's Silicone Cleaning Spray and rinsing with warm water for optimal results, then placing the device on the stand to dry. Never use cleaning products containing alcohol, petrol or acetone, as they may irritate the skin and damage the silicone.

FOREO MASKS

21. CAN I REUSE FOREO MASKS?

Each UFO™ Power Activated Mask and FOREO Sheet Mask is intended for single-use only. However, if one UFO™ 2 treatment wasn't enough, press the universal button again within 30 seconds after the end of the routine to start the same treatment over again!

22. DO I NEED TO WASH THE MASK ESSENCE OFF AFTER USING MY UFO™ 2 DEVICE?

No, Hyper-Infusion Technology is specifically engineered to help mask essence absorb deeper into the skin, faster. After the treatment, simply massage any remaining essence into your skin, or remove excess with a cotton pad.

23. HOW DO I KNOW IF MY FOREO MASK IS AUTHENTIC?

Using the FOREO For You app, you can scan the barcode on your mask packaging to verify its authenticity.

***NOTE:** Heat and cold affect different skincare products in different ways, and can compromise the ingredients in the formulas. (For optimal safety, we strongly recommend researching the effects of temperature on any non-FOREO mask before enhancing it with UFO™ 2.)

©2021 FOREO AB. ALL RIGHTS RESERVED.

EU DISTR.: FOREO AB, KARLAVÄGEN 41, 114 31 STOCKHOLM, SWEDEN.

US DISTR.: FOREO INC., 1525 E PAMA LANE, LAS VEGAS, NV 89119, USA.

AUS DISTR.: 6 PROHASKY STREET, PORT MELBOURNE, VIC, 3207 AUSTRALIA.

DESIGNED AND DEVELOPED BY FOREO SWEDEN.

MANUFACTURED FOR FOREO AB.

WWW.FOREO.COM