

FOREO  
SWEDEN



ESPADA™ 2

BLUE LED LIGHT SPOT TREATMENT



ONLINE MANUAL

# FULL USER MANUAL

---

## GETTING STARTED

Congratulations on taking the first step toward smarter skincare with ESPADA™ 2. Before you begin to enjoy all the benefits of sophisticated acne care in the comfort of your home, please take a few moments to look through this manual.

Please READ ALL INSTRUCTIONS BEFORE USE and utilize this product only for its intended use as described in this manual.

**⚠ WARNING:** NO MODIFICATION OF THIS EQUIPMENT IS ALLOWED.

## ESPADA™ 2 OVERVIEW

Featuring the dual-power of blue LED light & T-Sonic™ pulsations, ESPADA™ 2 works quickly & effectively to clear existing spots, blemishes & inflammation, and regulate sebum production to prevent future breakouts - leaving your skin looking clear, healthy & smooth.

## GETTING TO KNOW YOUR ESPADA™ 2

### FEATURES

#### BLUE LED LIGHT

415 nm wavelength penetrates the epidermis to kill spot-causing bacteria at the source.

#### SMART SKIN SENSOR

Ensures blue LED only activates when treatment window is in full contact with skin.

#### UNIVERSAL BUTTON

Turns the device on and off.

#### T-SONIC™ PULSATIONS

Massages the skin, stimulating circulation and helping to accelerate skin renewal.

#### BUILT-IN TIMER

Pulses to indicate when blemish has been treated, so you can move to the next spot.

#### TARGETING LIGHT

Helps to target the problem area, and aims blue LED at it with ultimate precision.

#### ANTIBACTERIAL SILICONE

Ultra-hygienic to reduce bacteria spreading from one area of skin to another.

#### CHARGING PORT

Up to 150 uses per full USB charge. 100% waterproof closure.

#### INDICATOR LIGHT

White indicates that device is on. Red indicates low battery.



# HOW TO USE YOUR ESPADA™ 2

## REGISTER YOUR WARRANTY

Before first use, download the FOREO For You app to register your device and activate the warranty by following these simple steps:

1. Download the FOREO app on your mobile phone by scanning the QR code.
2. Log in to your existing account, or sign up for a new one.
3. Add device (on the top of your screen).
4. Choose device series
5. Follow the instructions on the app and fill in the purchase information.

And your device is ready to use!



## HOW TO USE



### CLEANSE

Cleanse and dry your face and neck. For optimal results, we recommend using a FOREO LUNA™ Facial Cleansing Massager.



### APPLY

Apply FOREO ESPADA™ BHA+PHA Blemish Solution to each spot you wish to treat, if desired.



### TREAT

Press the universal button to turn on ESPADA™ 2. Aim it so the spot you wish to treat is in the center of the precise targeting light & hold on skin for 30 secs. You can treat the same spot up to 4 times in one treatment.

## CLEANING YOUR ESPADA™ 2

Always clean ESPADA™ 2 thoroughly after use. Clean the device with soap, rinse with warm water, and dry with a soft, lint-free cloth. After use, we recommend spraying the device with FOREO's Silicone Cleaning Spray and rinsing it with warm water for optimal results.

**NOTE:** Never use cleaning products containing alcohol, petrol, or acetone, as they may irritate the skin and damage the silicone.

## IMPORTANT FOR OPTIMUM SAFETY:

- If you have a severe skin condition or any medical concerns, please consult a dermatologist before use.
- Using ESPADA™ 2 should be comfortable – if you experience discomfort, discontinue use immediately and consult a physician.
- Exercise particular care when using ESPADA™ 2 around the eye area. Do not bring the device into contact with the eyelids or eyes themselves.
- For reasons of hygiene, we do not recommend sharing your ESPADA™ 2 with anyone else.
- Avoid leaving your ESPADA™ 2 in direct sunlight, and never expose it to extreme heat or boiling water.
- Close supervision is necessary when this device is used by, on, or near children, as well as those with reduced physical and mental abilities.
- Discontinue use if this product appears damaged in any way. This product contains no serviceable parts.
- Use this device only for its intended use as described in this manual. If you do not find the answer to your specific question, or if you have any other questions regarding the device's operation, please visit [www.foreo.com](http://www.foreo.com).

## TROUBLESHOOTING

*Precautions to be taken in the event of changes in the performance of ESPADA™ 2.*

### **If ESPADA™ 2 is not activated when pressing the universal button:**

- Battery is empty. Charge your ESPADA™ 2 via the provided USB charging cable for 1.5 hrs until fully charged. A full charge will last up to 150 thirty-second treatments.

### **If ESPADA™ 2 cannot be switched off and/or interface buttons do not respond:**

- Microprocessor is temporarily malfunctioning. Press and hold the universal button to restart the device. If ESPADA™ 2's indicator light does not blink when the USB charging cable is connected:
- Battery is fully charged and ready for up to 150 treatments.
- Battery is drained and will require several minutes to acknowledge USB charging cable.
- USB charging cable is not properly connected. Check electrical outlet / USB socket and charging port.

## WARRANTY TERMS & CONDITIONS

### **REGISTER WARRANTY**

To activate your 2-Year Limited Warranty, register through the FOREO For You app, or visit [foreo.com/product-registration](https://foreo.com/product-registration) for more information.

### **2-YEAR LIMITED WARRANTY**

FOREO warrants this device for a period of TWO (2) YEARS (except in countries where national law requires a longer minimum guarantee) after the original date of purchase against defects due to faulty workmanship or materials arising from Normal Use of the device. The warranty covers working parts that affect the function of the device. It does NOT cover cosmetic deterioration caused by fair wear and tear, or damage caused by accident, misuse or neglect. Any attempt to open or take apart the device (or its accessories) will void the warranty.

If you discover a defect and notify FOREO during the warranty period, FOREO will, at its discretion, replace the device free of charge. Claims under warranty must be supported by reasonable evidence that the date of the claim is within the warranty period. To validate your warranty, please keep your original purchase receipt together with these warranty conditions for the duration of the warranty period.

To claim your warranty, you must log in to your account at [www.foreo.com](https://www.foreo.com) and then select the option to make a warranty claim. Shipping costs are non-refundable. This undertaking is in addition to your statutory rights as a consumer and does not affect those rights in any way.

## DISPOSAL INFORMATION

Disposal of old electronic equipment (applicable in the EU and other European countries with separate waste collection systems).




The crossed-out dust bin symbol indicates that this device should not be treated as household waste, but rather be brought to the appropriate collection point for recycling of electrical and electronic equipment. By ensuring this device is disposed of correctly, you will help prevent the potential negative consequences for the environment and human health which could be caused by inappropriate waste handling of the product. The recycling of materials will also help conserve natural resources.

For more information about the recycling of your device, please contact your local household waste disposal service or your place of purchase.

# BATTERY REMOVAL

**NOTE:** This process is not reversible. Opening the device will void its warranty. This action must only be undertaken when the device is ready to be disposed of.

Because this device contains a lithium-ion battery, the battery must be removed before disposal and should not be thrown away with household waste. To remove the battery, open the inner plastic shell after removing silicone outer layer and remove the battery to be disposed of in accordance with your local environmental regulations. Wear gloves during this process for your safety. Detailed visual instructions are provided below.



# SPECIFICATIONS

<b>MATERIALS:</b>	Body-safe silicone, PC + ABS	<b>USAGE:</b>	Up to 150 uses per charge
<b>COLOR:</b>	Cobalt Blue / Fuchsia / Pearl Pink	<b>STANDBY:</b>	180 days
<b>SIZE:</b>	H140 x W36 x D31 mm	<b>FREQUENCY:</b>	103 Hz
<b>WEIGHT:</b>	87 g	<b>MAX NOISE LEVEL:</b>	<50 dB
<b>BATTERY:</b>	420 mAh	<b>INTERFACE:</b>	1-button
<b>TEMPERATURE:</b>	5 - 45 degrees celsius		

**DISCLAIMER:** Users of this device do so at their own risk. Neither FOREO nor its retailers assume any responsibility or liability for any injuries or damages, physical or otherwise, resulting directly or indirectly from the use of this device. Furthermore, FOREO reserves the right to revise this publication and to make changes from time to time in the contents thereof without obligation to notify any person of such revision or changes.

**CAUTION:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Model may be changed for improvements without notice.

©2023 FOREO AB. ALL RIGHTS RESERVED. EU DISTR.: FOREO AB, KARLAVÄGEN 41, 114 31 STOCKHOLM, SWEDEN. US DISTR.: FOREO INC., 1525 E PAMA LANE, LAS VEGAS, NV 89119, USA. AUS DISTR.: 6 PROHASKY STREET, PORT MELBOURNE, VIC, 3207 AUSTRALIA. DESIGNED AND DEVELOPED BY FOREO SWEDEN. MANUFACTURED FOR FOREO AB.